

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY  
DEVICE ONLY TEMPLATE**

**A. 510(k) Number:**

k042243

**B. Purpose for Submission:**

New Device

**C. Analyte:**

Albumin

**D. Type of Test:**

Quantitative

**E. Applicant:**

Diagnostic Chemicals Limited

**F. Proprietary and Established Names:**

Microalbumin Multi-Calibrator Set

**G. Regulatory Information:**

1. Regulation section:  
21 CFR 862.1150 Calibrator
2. Classification:  
II
3. Product Code:  
JIX
4. Panel:  
75

**H. Intended Use:**

1. Intended use(s):  
For IN VITRO diagnostic use as a calibrator for the DCL Microalbumin Assay for quantitation of albumin in human urine.
2. Indication(s) for use:  
For IN VITRO diagnostic use as a calibrator for the DCL Microalbumin Assay for quantitation of albumin in human urine.
3. Special condition for use statement(s):  
For Prescription Use Only

4. Special instrument Requirements:  
N/A

**I. Device Description:**

The Microalbumin Multi-Calibrator set contains 6 levels of 1 mL of liquid solutions corresponding to the following values: 1, 5, 10, 50, 100 and 300 mg/L. The human donor unit used to produce this control was tested by FDA accepted methods and found free of Hepatitis B Surface Antigen (HBsAG), antibody to Hepatitis C (HCV) and antibody to HIV.

**J. Substantial Equivalence Information:**

1. Predicate device name(s):  
Kamiya Biomedical Company Microalbumin Calibrator Set
2. Predicate K number(s):  
K991166
3. Comparison with predicate:

<b>Similarities</b>		
<b>Item</b>	<b>Device</b>	<b>Predicate</b>
Intended Use	For IN VITRO diagnostic use as a calibrator for the DCL Microalbumin Assay for quantitation of albumin in human urine.	The K-Assay microalbumin Calibrator Set is used for the calibration of the K-Assay Microalbumin immunoturbimetric assay for quantitating albumin in urine specimens.
Sample	Urine	Urine
Matrix	Liquid	Liquid

**K. Standard/Guidance Document Referenced (if applicable):**

ICH Biocompatibility Guidance

**L. Test Principle:**

Not Applicable

**M. Performance Characteristics (if/when applicable):**

1. Analytical performance:
  - a. *Precision/Reproducibility:*  
Not Applicable
  - b. *Linearity/assay reportable range:*  
Not Applicable
  - c. *Traceability (controls, calibrators, or method):*  
Scripps Lab (Cat # A0213) Human Serum

Stability

Once opened, the calibrators are stable for 30 days when stored at 2-8 °C. An accelerated study was conducted on the Microalbumin Calibrator set to support the shelf life of 12 months. An on-going real time study is being conducted and the Microalbumin Calibrator set is checked at its half life and at 13 months.

#### Value Assignment

An inhouse standard is used to prepare a calibration curve by diluting it to the desired concentrations. The values of the calibrator concentrations are assigned from multiple readings using DCL Microalbumin reagent on a Hitachi 717 or equivalent analyzer. The acceptable CV % for each sample is 3 percent. The recovery rate of each concentration of the diluted primary standard is 1.00 +/- 0.03. The calibrator lots are validated using the internal standard, serum with a known value, and the previous lot of calibrators. The new lot is acceptable if the sample results fall within 100 +/- 5% of their known values.

d. *Detection limit:*

N/A

e. *Analytical specificity:*

N/A

f. *Assay cut-off:*

N/A

#### 2. Comparison studies:

a. *Method comparison with predicate device:*

N/A

b. *Matrix comparison:*

N/A

#### 3. Clinical studies:

a. *Clinical sensitivity:*

N/A

b. *Clinical specificity:*

N/A

c. *Other clinical supportive data (when a and b are not applicable):*

N/A

#### 4. Clinical cut-off:

N/A

#### 5. Expected values/Reference range:

The concentrations of the calibrators in the Microalbumin Multi-Calibrator Set have lot specific values that range from 0 to 300 mg/L.

**N. Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.