

## **510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY**

**A. 510(k) Number:**

k061970

**B. Purpose for Submission:**

Modification to the manufacturer's existing device including changes to a strictly mono-biotinylated monoclonal antibody (the antibody itself remains unchanged), deletion of an application for the Elecsys 1010 analyzer and modifications to specifications for linearity, intra-assay and inter-assay precision.

**C. Measurand:**

Quantitative determination of IgE in human serum and plasma

**D. Type of Test:**

Immunoassay

**E. Applicant:**

Roche Diagnostics.

**F. Proprietary and Established Names:**

Elecsys IgE II immunoassay

**G. Regulatory Information:**

1. Regulation section:  
21 CFR 866.5510, Immunoglobulins A, G, M, D and E Immunological test system
2. Classification:  
Class II
3. Product code:  
JHR, Radioimmunoassay, Immunoglobulins (D, E)
4. Panel:  
Immunology (82)

**H. Intended Use:**

1. Intended use(s):  
Immunoassay for the in vitro quantitative determination of immunoglobulin E in human serum and plasma. Determination of total IgE is useful as an aid in the diagnosis of allergic diseases.
2. Indication(s) for use:  
Same as Intended Use
3. Special conditions for use statement(s):  
For Prescription use only
4. Special instrument requirements:

**I. Device Description:**

The Elecsys IgE II immunoassay is a two step sandwich immunoassay with streptavidin microparticles and electrochemiluminescence detection. Results are determined using a calibration curve that is generated specifically on each instrument by a 2 point calibration and a master curve provided with the reagent bar code.

**J. Substantial Equivalence Information:**

1. Predicate device name(s):

- Elecsys IgE immunoassay
2. Predicate 510(k) number(s):  
(K984326, K961481/A003)
  3. Comparison with predicate:

Topic	Elecsys IgE (K984326, K961481/A003)	Elecsys IgE II (Modified Device)
Intended use	Immunoassay for the in vitro quantitative determination of immunoglobulin E in human serum and plasma. Determination of total IgE is useful as an aid in the diagnosis of allergic diseases.	Same
Analyzers	Roche Elecsys 1010/2010 and MODULAR ANALYTICS E170 (Elecsys module) immunoassay analyzers.	Roche Elecsys 2010 and MODULAR ANALYTICS E170 (Elecsys module) immunoassay analyzers.
Traceability	Assay standardized against the 2 <sup>nd</sup> IRP WHO Reference Standard 75/502	Same
Assay Protocol	Sandwich	Same
Sample Type	Serum and plasma	Same
Calibrator, Calibration Verification & Controls	IgE CalSet IgE CalCheck PreciControl Universal	Same
Measuring Range	0.100 – 2500 IU/ml	Same
Analytical Sensitivity	0.10 IU/ml	Same
Functional Sensitivity	0.50 IU/ml	Same
Composition	R1: 2.4 mg/L AB-Bi, buffer, preservative R2: 4.8 mg/L AB-Ru, buffer, preservative M: 0.72 mg/ml streptavidin-coated microparticles, preservative	R1: 2.5 mg/L (mono-Bi)AB-Bi, buffer, preservative R2: 5.5 mg/L AB-Ru, buffer, preservative M: 0.72 mg/ml streptavidin-coated microparticles, preservative
Dilution Recommendation	Concentration of diluted samples must be > 60 IU/mL	Concentration of diluted samples must be > 125 IU/mL

**K. Standard/Guidance Document Referenced (if applicable):**

**L. Test Principle:**

Same

**M. Performance Characteristics (if/when applicable):**

1. Analytical performance:
  - a. Precision/Reproducibility:

**Special 510(k): Device Modification - Elecsys IgE II**  
**Immunoassay, continued**

Topic	Elecsys IgE (K984326, K961481/A003)	Elecsys IgE II (Modified Device)
Interferences	No affect up to, bilirubin < 37 mg/dl hemoglobin < 1.1 g/dl triglycerides < 2200 mg/dl biotin < 100 ng/ml	No affect up to, Bilirubin - Same hemoglobin < 0.1 g/dl triglycerides - Same biotin - Same
Precision	<i>Elecsys 1010 / 2010:</i> Intra-assay HS1 3.6% CV @ 5.18 IU/mL HS2 3.2% CV @ 398 IU/mL HS3 2.4% CV @ 1010 IU/mL Total: HS1 4.2% CV @ 5.18 IU/mL HS2 3.9% CV @ 398 IU/mL HS3 3.1% CV @ 1010 IU/mL	<i>Elecsys 2010:</i> Intra-assay: HS1 4.1% CV @ 32.7 IU/mL HS2 2.4% CV @ 265 IU/mL HS3 2.6% CV @ 1295 IU/mL Total: HS1 5.1% CV @ 32.7 IU/mL HS2 3.8% CV @ 265 IU/mL HS3 3.9% CV @ 1295 IU/mL

b. *Linearity/assay reportable range:*

Not applicable

c. *Traceability, Stability, Expected values (controls, calibrators, or methods):*

d. *Detection limit:*

Same

e. *Analytical specificity:*

Same

f. *Assay cut-off:*

Not applicable.

2. Comparison studies:

a. *Method comparison with predicate device:*

b. *Matrix comparison:*

Not applicable

3. Clinical studies:

a. *Clinical Sensitivity:*

No change.

b. *Clinical specificity:*

No change.

c. *Other clinical supportive data (when a. and b. are not applicable):*

Not applicable.

4. Clinical cut-off:

Not applicable.

5. Expected values/Reference range:

N. Software:

**O. Proposed Labeling:**

**P. Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.