

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
DEVICE ONLY TEMPLATE**

A. 510(k) Number:

#K033491

B. Purpose for Submission:

New device

C. Analyte:

D-Dimer

D. Type of Test:

Quantitative

E. Applicant:

Roche Diagnostics Corporation

F. Proprietary and Established Names:

CARDIAC D-Dimer Assay; Fibrin Split Products, Antigen, antiserum, Control

G. Regulatory Information:

1. Regulation section:
CFR Section 864.7320 - Fibrin Degradation Products Assay
2. Classification:
Class II
3. Product Code:
DAP; GHH
4. Panel:
Hematology (81)

H. Intended Use:

1. Intended use(s):
The Roche CARDIAC D-Dimer Assay is intended for the quantitative determination of D-Dimer in anticoagulated venous whole blood with the CARDIAC reader instrument.
2. Indication(s) for use:
Same as the Intended Use.
3. Special condition for use statement(s):
N/A
4. Special instrument Requirements:
Roche Diagnostics Corp. CARDIAC Reader (#K000784)

I. Device Description:

The CARDIAC D-Dimer Assay consists of a test strip in plastic casing, with a well for sample application; viewing window at the detection zone; a lot-specific coding chip for each batch of test strips; and Control, Levels I/II. Testing is initiated by the addition of whole blood to the test strip, where red cells are separated from the plasma. A monoclonal antibody (MAB) sandwich complex is formed and bound to a streptavidin band in the detection zone, where the gold-labelled antibody forms a reddish-purple color signal. That signal is determined by the CARDIAC reader.

J. Substantial Equivalence Information:

1. Predicate device name(s):
Roche Diagnostics Corp. Tina-quant® D-Dimer Test System
2. Predicate K number(s):
#K002706
3. Comparison with predicate:

Similarities		
Item	Device	Predicate
Sample anticoagulant	Heparin	Same
Assay type	Quantitative immunoassay that utilizes mouse MAB's	Same
Reagents	Includes buffer	Same
Differences		
Item	Device	Predicate
Sample type	Heparinized plasma separated from venous whole blood during the test procedure	Pre-prepared citrated or heparinized plasma
Principle	Sandwich assay with color indicator	Turbidometric latex assay
Reagents	Biotinylated and gold-labelled anti-D-Dimer antibody	Latex particles coated with anti-human D-Dimer MAB
Measuring range	0.1 – 4.0 ug/mL	0.15 – 9.0 ug/mL

K. Standard/Guidance Document Referenced (if applicable):

N/A

L. Test Principle:

The Roche CARDIAC D-Dimer Assay is an immunoassay based upon the dual monoclonal antibody (MAB) sandwich principle. It uses a poly-streptavidin-biotin capture system with a gold sol particle label. D-Dimer in the sample combines with the antibodies conjugated to the gold sol particles to form a “sandwich”. This, in turn, combines with the poly-streptavidin-biotin, which is immobilized in a line across the result window. The speed and intensity of color formation is related to the concentration of D-Dimer in the sample; and results are determined on the CARDIAC Reader instrument.

M. Performance Characteristics (if/when applicable):1. Analytical performance:a. *Precision/Reproducibility:*

(2) to (3) lots of Levels I/II quality control were run at each site.

Within-run – performed at (4) sites and run X (10) on each analyzer. Level I %CV's ranged 7.8 – 14 %. Level II %CV's ranged 8.0 – 18%.

Day-to-day – performed at (2) sites, over 20 - 36 days on each analyzer. Level I %CV's ranged 12.3 – 14%. Level II %CV's ranged 10.5 – 17.1%.

Acceptance criteria: $\leq 20\%$

b. *Linearity/assay reportable range:*

0.1 – 4 mg/L, determined from method comparison data.

c. *Traceability (controls, calibrators, or method):*

Calibrated to the Roche Tina-quant D-Dimer method by using a primal master lot to establish a calibration curve.

d. *Detection limit:*

0.1 mg/L

e. *Analytical specificity:*

Various levels of D-Dimer tested for endogenous substances (hematocrit, hemoglobin, bilirubin, lipids, and HAMA); and exogenous substances (biotin, heparin, Vitamin C, aspirin and prescription drugs). Fibrinogen and fragments D, E, X and Y were also tested. Rheumatoid factor was not tested.

f. *Assay cut-off:*

N/A

2. Comparison studies:

a. *Method comparison with predicate device:*

Performed at (6) clinical sites, including (5) point-of-care (POC) sites. Studies were done on patients suspected of thromboembolic venous diseases, peri- or post-OR patients and healthy donors. Three lots of CARDIAC test strips were run against the Tina-quant assay. The following regression data were obtained:

Lot 1, (6) sites - (N = 464); $y = -0.04 + 1.09x$; $r = 0.89$

Lot 2, (6) sites - (N = 463); $y = -0.02 + 1.08x$; $r = 0.90$

Lot 3, (4) sites - (N = 92); $y = -0.02 + 1.30x$; $r = 0.90$

b. *Matrix comparison:*

Ammonium and lithium heparins were tested on samples (N = 164) with the CARDIAC and Tina-quant D-Dimer assays. Regression data: $y = 0.02 + 1.118x$; $r = 0.908$.

3. Clinical studies:

a. *Clinical sensitivity:*

N/A

b. *Clinical specificity:*

N/A

c. *Other clinical supportive data (when a and b are not applicable):*

N/A

4. Clinical cut-off:

N/A

5. Expected values/Reference range:

A normal range study was performed on normal individuals (N = 130), consisting of 93 females and 37 males, ranging in age from 19 - 94 years. The mean value was 0.50 mg/L.

N. Conclusion:

The submitted material in this premarket notification is complete and supports a substantial equivalence decision.