

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY  
ASSAY ONLY TEMPLATE**

**A. 510(k) Number:**

K062203

**B. Purpose for Submission:**

Modification to the intended use and reagent formulation (a pH change of a buffer)

**C. Measurand:**

D-dimer

**D. Type of Test:**

Immunoturbidimetric

**E. Applicant:**

Roche Diagnostics

**F. Proprietary and Established Names:**

Tina-Quant D-dimer

**G. Regulatory Information:**

1. Regulation section:

21 CFR 864.7320

2. Classification:

Class II

3. Product code:

GHH

4. Panel:

81 Hematology

## **H. Intended Use:**

### 1. Intended use(s):

The Roche diagnostics Tina-Quant D-dimer assay is an immunoturbidimetric assay for the in vitro quantitative determination of fibrin degradation products including d-dimer and x-oligomers in human plasma. In conjunction with a non-high clinical probability assessment, a normal ( $<0.5 \mu\text{g FEU/mL}$ ) result excludes deep vein thrombosis (DVT) and pulmonary embolism (PE) with high sensitivity.

### 2. Indication(s) for use:

### 3. Special conditions for use statement(s):

### 4. Special instrument requirements:

Roche automated clinical chemistry analyzers

## **I. Device Description:**

The Roche Tina-quant D-dimer is an immunoturbidimetric assay consisting of a ready to use tris buffer and a ready to use anti-d-dimer latex suspension.

## **J. Substantial Equivalence Information:**

### 1. Predicate device name(s):

a. Roche Diagnostics Tina-Quant D-Dimer Test System

b. Biomerieux Vidas D-Dimer Exclusion Assay

### 2. Predicate 510(k) number(s):

a. K030740

b. K040822

### 3. Comparison with predicate:

Similarities			
Item	Device	Predicate (a)	Predicate (b)
Intended Use	Quantitative d-dimer and in conjunction with PTP, exclude DVT and PE		same
Quality control	Bi-level	Same	same

Item	Device	Predicate (a)	Predicate (b)
Intended Use	Quantitative d-dimer and in conjunction with PTP, exclude DVT and PE	Quantitative d-dimer	
Technology	Immunoturbidimetric	Same	Enzyme Linked Fluorescent Assay
Sample Matrix	Citrated Plasma Lithium Heparin Plasma	Same	Citrated Plasma
Buffer	Anti-D-dimer latex suspension in Tris Buffer pH7.2	Anti-D-Dimer latex suspension in Tris buffer pH 8.2	NA

**K. Standard/Guidance Document Referenced (if applicable):**

**L. Test Principle:**

The Tina-quant D-Dimer Test System is based on a particle enhanced immunoturbidimetric method. Latex particles of uniform size are coated with monoclonal antibodies (F(ab')<sub>2</sub> fragments) to the D-dimer epitope. The antigen/antibody complexes produced by the addition of samples containing D-dimer lead to an increase in the turbidity of the test reactants, which can be determined

turbidmetrically.

**M. Performance Characteristics (if/when applicable):**

1. Analytical performance:

a. *Precision/Reproducibility:*

<u>Within-run</u>	Roche/Hitachi
0.19µg FEU/mL	7.3%
0.86 µg FEU/mL	1.7%
5.11 FEU/mL	0.8%
	COBAS Integra
0.279 µg FEU/mL	6.9%
2.88 µg FEU/mL	1.1%

b. *Linearity/assay reportable range:*

c. *Traceability, Stability, Expected values (controls, calibrators, or methods):*

d. *Detection limit:*

Roche/Hitachi:	COBAS Integra:
0.04 µg FEU/mL	<0.08 µg FEU/mL

e. *Analytical specificity:*

f. *Assay cut-off:*

2. Comparison studies: N/A

a. *Method comparison with predicate device:*

b. *Matrix comparison:*

3. Clinical studies:

a. *Clinical Sensitivity:*

DVT Exclusion

Tina-Quant® D-dimer was used in a multicenter management study involving 812 outpatients with suspected DVT. Patient pre-test probability (PTP) was assessed following Wells probability score. Patients were classified as high ( $\geq 3$ ) or non-high ( $< 3$ ), and the Tina-Quant® D-dimer test was then performed using a cutoff of 0.5  $\mu\text{g}$  FEU/mL. Those patients with a normal (negative) d-dimer, and a non-high pretest probability had no further diagnostic testing and were followed for 3 months for development of DVT. 1 of 176 such patients developed DVT during the follow-up period. Sensitivity, NPV and Failure rates were determined as follows:

Sensitivity:	99.3%
NPV	99.4%
Failure Rate:	0.6%

PE Exclusion

Tina-Quant® D-Dimer was used in a management study involving 168 outpatients with suspected PE. Patient pre-test probability (PTP) was assessed following Wells probability score for PE. Patients were classified as having a low, moderate, or high pretest probability of PE, and the Tina-Quant® D-dimer test was then performed using a cutoff of 0.5  $\mu\text{g}$  FEU/mL. Those patients with a normal (negative) d-dimer, and a non-high (low or moderate) pretest probability had no further diagnostic testing and were followed for 3 months for development of PE. No patients developed PE during the follow-up period. Sensitivity, NPV and Failure rates were determined as follows:

Sensitivity:	100%
NPV	100%
Failure Rate:	0%

b. *Clinical specificity:*

c. Other clinical supportive data (when a. and b. are not applicable):

4. Clinical cut-off:

N/A

5. Expected values/Reference range:

N/A

**N. Proposed Labeling:**

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

**O. Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.